



- en** Addendum **da** Tilføjelse **de** Addendum **es** Apéndice **fi** Lisäys
- fr** Addendum **it** Addendum **nl** Aanvulling **no** Tillegg **pt** Adenda
- sv** Tillägg

3.2 Options / Ekstraustyr / Optionen / Opciones / Lisävarusteet / Options / Opzioni / Opties / Ekstraustyr / Opções / Tillval

- Extension loops / Forlængerløkker / Verlängerungsschlaufen / Anillas de extensión / Jatkesilmukat / Boucles d'extension / Occhielli di estensione / Verlenglusen / Forlængelsesløkker / Presilhas de extensão / Förlängningsöglor
- Loop marker/ Loop-markør / Schlaufenmarkierung / Marcador de anillas / Silmukkamerkki / Marqueur de boucle / Indicatore di occhielli / Lusmarkering / Løkkemarkør / Marcador de presilha / Ögelmarkör

7 Technical Data / Tekniske Data / Technische Daten / Datos Técnicos / Tekniset Tiedot / Caractéristiques Techniques / Dati Tecnici / Technische Specificaties / Tekniske Data / Características Técnicas / Teknisk Data

Dimensions / Mål / Abmessung / Dimensiones / Mitat / Dimenszioni / Afmetingen / Dimensjoner / Dimensões / Mått

Variant / Ausführung / Variante / Versio	[mm]	XS	S	M	L	XL	
Stand Assist	A	260	260	330	360	360	
	B	840	840	895	950	1020	
	C	900	950	1170	1270	1810	
Transfer Stand Assist	A	N/A	260	330	360	360	
	B	N/A	840	895	950	1020	
	C	N/A	950	1170	1270	1810	
	D	N/A	965	1060	1160	1260	

60147943-A 2024-09-25 1/1



Invacare UK Operations Limited
Unit 4, Pencoed Technology Park,
Pencoed
Bridgend CF35 5AQ
UK



Invacare Portugal, Lda
Rua Estrada Velha 949
4465-784 Leça do Balio
Portugal



Yes, you can.®